

Conjugated CBL (Phospho-Tyr774) Rabbit mAb

Catalog No: #C14270

Package Size: #C14270-AF350 100ul #C14270-AF405 100ul #C14270-AF488 100ul #C14270-AF555 100ul #C14270-AF594 100ul
#C14270-AF647 100ul #C14270-AF680 100ul #C14270-AF750 100ul #C14270-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated CBL (Phospho-Tyr774) Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	Phospho-CBL (Y774) Antibody detects endogenous levels of total Phospho-CBL (Y774)
Immunogen Description	A synthesized peptide derived from human Phospho-CBL (Y774)
Conjugates	AF Dye
Modification	Phospho
Other Names	Casitas B lineage lymphoma proto oncogene; CBL2; E3 ubiquitin protein ligase CBL; Oncogene CBL2; Proto oncogene c CBL; RING finger protein 55; RNF55; Signal transduction protein CBL;
Accession No.	Uniprot:P22681
Calculated MW	Predicted band size: 100 kDa
SDS-PAGE MW	Observed band size: 120 kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution::

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.