

Conjugated TAOK1/2/3 (Phospho-Ser181/Ser181/Ser177) Rabbit mAb

Catalog No: #C14277

Orders: order@signalwayantibody.com

Package Size: #C14277-AF350 100ul #C14277-AF405 100ul #C14277-AF488 100ul #C14277-AF555 100ul #C14277-AF594 100ul

#C14277-AF647 100ul #C14277-AF680 100ul #C14277-AF750 100ul #C14277-Biotin 100ul

Description

Product Name	Conjugated TAOK1/2/3 (Phospho-Ser181/Ser181/Ser177) Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Phospho-TAOK1/2/3 (S181 + S181 + S177) Antibody detects endogenous levels of total Phospho-TAOK1/2/3 (S181 + S181 + S177)
Immunogen Description	A synthesized peptide derived from human Phospho-TAOK1/2/3 (S181 + S181 + S177)
Conjugates	AF Dye
Modification	Phospho
Other Names	TAOK1; TAOK2; TAOK3;
Accession No.	Uniprot:Q7L7X3/Q9UL54/Q9H2K8
Calculated MW	Predicted band size: 116 kDa
SDS-PAGE MW	110-140 kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4ℳ Cin dark for 6 months

Application Details

Suggested Dilution;;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.