

Conjugated p130(Phospho-Thr986) Rabbit mAb

Catalog No: #C14282

Package Size: #C14282-AF350 100ul #C14282-AF405 100ul #C14282-AF488 100ul #C14282-AF555 100ul #C14282-AF594 100ul
 #C14282-AF647 100ul #C14282-AF680 100ul #C14282-AF750 100ul #C14282-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	Conjugated p130(Phospho-Thr986) Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse
Specificity	Phospho-p130(T986) Antibody detects endogenous levels of total Phospho-p130(T986)
Immunogen Description	A synthesized peptide derived from human Phospho-p130(T986)
Conjugates	AF Dye
Modification	Phospho
Other Names	130 kDa retinoblastoma associated protein; IB4; p130; PRB2; Rb2; Rbl2; RBR2; Retinoblastoma like 2; Retinoblastoma-related protein 2; RIC128;
Accession No.	Uniprot:Q08999
Calculated MW	128kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution;;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.