

Conjugated c-Jun/JunD(Phospho-Ser73+Ser100) Rabbit mAb

Catalog No: #C14284

Orders: order@signalwayantibody.com

Package Size: #C14284-AF350 100ul #C14284-AF405 100ul #C14284-AF488 100ul #C14284-AF555 100ul #C14284-AF594 100ul

#C14284-AF647 100ul #C14284-AF680 100ul #C14284-AF750 100ul #C14284-Biotin 100ul

Description

Product Name	Conjugated c-Jun/JunD(Phospho-Ser73+Ser100) Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse
Specificity	Phospho-c-Jun (S73)+JunD (S100) Antibody detects endogenous levels of total Phospho-c-Jun (S73)+JunD (S100)
Immunogen Description	A synthesized peptide derived from human Phospho-c-Jun (S73)+JunD (S100)
Conjugates	AF Dye
Modification	Phospho
Other Names	JunD; c-Jun;
Accession No.	Uniprot:P05412,P17535
Calculated MW	Predicted band size: 36 kDa
SDS-PAGE MW	38-40 kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4 Λ Cin dark for 6 months

Application Details

Suggested Dilution;;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.