

MCU(Phospho-Ser92) Conjugated Antibody

Catalog No: #C14297

Package Size: #C14297-AF350 100ul #C14297-AF405 100ul #C14297-AF488 100ul #C14297-AF555 100ul #C14297-AF594 100ul
 #C14297-AF647 100ul #C14297-AF680 100ul #C14297-AF750 100ul #C14297-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MCU(Phospho-Ser92) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Phospho-MCU(S92) Antibody detects endogenous levels of total Phospho-MCU(S92)
Conjugates	AF Dye
Other Names	C10orf42; Calcium uniporter protein mitochondrial; Ccdc109a; HsMCU; Mitochondrial Calcium Uniporter;
Accession No.	Q8NE86
Calculated MW	30kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

Mitochondrial inner membrane calcium uniporter that mediates calcium uptake into mitochondria (PubMed:21685888, PubMed:21685886, PubMed:23101630, PubMed:22904319, PubMed:23178883, PubMed:22829870, PubMed:22822213, PubMed:24332854, PubMed:23755363, PubMed:26341627).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.