

# Conjugated PTBP2 Rabbit mAb

Catalog No: #C56010

Package Size: #C56010-AF350 100ul #C56010-AF405 100ul #C56010-AF488 100ul #C56010-AF555 100ul #C56010-AF594 100ul  
#C56010-AF647 100ul #C56010-AF680 100ul #C56010-AF750 100ul #C56010-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                               |
|-----------------------|---------------------------------------------------------------------------------------------------------------|
| Product Name          | Conjugated PTBP2 Rabbit mAb                                                                                   |
| Host Species          | Recombinant Rabbit                                                                                            |
| Clonality             | Monoclonal                                                                                                    |
| Purification          | Affinity-chromatography                                                                                       |
| Species Reactivity    | Human Mouse Rat                                                                                               |
| Specificity           | PTBP2 Antibody detects endogenous levels of total PTBP2                                                       |
| Immunogen Description | A synthesized peptide derived from human PTBP2                                                                |
| Conjugates            | AF Dye                                                                                                        |
| Other Names           | brPTB; MIBP; nPTB; Polypyrimidine tract binding protein 2; PTB like protein; PTBP 2Ptbp2; Splicing regulator; |
| Accession No.         | Uniprot:Q9UKA9                                                                                                |
| Calculated MW         | 57kDa                                                                                                         |
| Formulation           | Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.                                                    |
| Storage               | Store at 4A Cin dark for 6 months                                                                             |

## Application Details

Suggested Dilution.;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.