

Conjugated GM130 Rabbit mAb

Catalog No: #C56019

Package Size: #C56019-AF350 100ul #C56019-AF405 100ul #C56019-AF488 100ul #C56019-AF555 100ul #C56019-AF594 100ul
#C56019-AF647 100ul #C56019-AF680 100ul #C56019-AF750 100ul #C56019-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated GM130 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat Monkey Cow Dog
Specificity	GM130 Antibody detects endogenous levels of total GM130
Immunogen Description	A synthesized peptide derived from human GM130
Conjugates	AF Dye
Other Names	GM130; Gm130 autoantigen; GOLGA 2; Golga2; Golgi autoantigen; Golgin 95;
Accession No.	Uniprot:Q08379
Calculated MW	Predicted band size: 113 kDa
SDS-PAGE MW	Observed band size: 130 kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution:;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.