

# Conjugated MHC class 1 Rabbit mAb

Catalog No: #C56034

Package Size: #C56034-AF350 100ul #C56034-AF405 100ul #C56034-AF488 100ul #C56034-AF555 100ul #C56034-AF594 100ul #C56034-AF647 100ul #C56034-AF680 100ul #C56034-AF750 100ul #C56034-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                                       |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Conjugated MHC class 1 Rabbit mAb                                                                                                     |
| Host Species          | Recombinant Rabbit                                                                                                                    |
| Clonality             | Monoclonal                                                                                                                            |
| Purification          | Affinity-chromatography                                                                                                               |
| Species Reactivity    | Human                                                                                                                                 |
| Specificity           | MHC class 1 Antibody detects endogenous levels of total MHC class 1                                                                   |
| Immunogen Description | A synthesized peptide derived from human MHC class 1                                                                                  |
| Conjugates            | AF Dye                                                                                                                                |
| Other Names           | Aw-68;HLA class I histocompatibility antigen, A-28 alpha chain;MHC class I antigen A*68;HLA-A; MHC class I antigen HLA A heavy chain; |
| Accession No.         | Uniprot:P01891                                                                                                                        |
| Calculated MW         | 41kDa                                                                                                                                 |
| Formulation           | Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.                                                                            |
| Storage               | Store at 4ℳ Cin dark for 6 months                                                                                                     |

## Application Details

Suggested Dilution:;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.