

Conjugated DAP Kinase 1 Rabbit mAb

Catalog No: #C56044

Package Size: #C56044-AF350 100ul #C56044-AF405 100ul #C56044-AF488 100ul #C56044-AF555 100ul #C56044-AF594 100ul
 #C56044-AF647 100ul #C56044-AF680 100ul #C56044-AF750 100ul #C56044-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated DAP Kinase 1 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human;Mouse;Rat
Specificity	DAP Kinase 1 Antibody detects endogenous levels of total DAP Kinase 1
Immunogen Description	A synthesized peptide derived from human DAP Kinase 1
Conjugates	AF Dye
Other Names	DAK1; DAP K1; DAP kinase 1; DAPK 1; DAPK; DAPK1; DAPK1_HUMAN; Death Associated Protein Kinase 1; Death-associated protein kinase 1;
Accession No.	Uniprot:P53355
Calculated MW	160 kDa
SDS-PAGE MW	160 kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4ℳ Cin dark for 6 months

Application Details

Suggested Dilution;;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.