

Conjugated RAB11A Rabbit mAb

Catalog No: #C56050

Package Size: #C56050-AF350 100ul #C56050-AF405 100ul #C56050-AF488 100ul #C56050-AF555 100ul #C56050-AF594 100ul
#C56050-AF647 100ul #C56050-AF680 100ul #C56050-AF750 100ul #C56050-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated RAB11A Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	RAB11A Antibody detects endogenous levels of total RAB11A
Immunogen Description	A synthesized peptide derived from human RAB11A
Conjugates	AF Dye
Other Names	MGC1490; Rab 11; Rab 11A; Rab-11; RAB11 A; RAB11; Rab11a; RAB11A member RAS oncogene family; Ras related protein Rab 11A; Ras related protein Rab11A; Ras-related protein Rab-11A; RB11A; YL 8; YL8;
Accession No.	Uniprot:P62491
Calculated MW	24kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution;;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.