

Conjugated CaMKII Rabbit mAb

Catalog No: #C56058

Package Size: #C56058-AF350 100ul #C56058-AF405 100ul #C56058-AF488 100ul #C56058-AF555 100ul #C56058-AF594 100ul
#C56058-AF647 100ul #C56058-AF680 100ul #C56058-AF750 100ul #C56058-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated CaMKII Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	CaMKII Antibody detects endogenous levels of total CaMKII
Immunogen Description	A synthesized peptide derived from human CaMKII
Conjugates	AF Dye
Other Names	CAMK2; Camk2a; CAMK2B; CAMK2D; CAMK2G; CAMKA;
Accession No.	Uniprot:Q9UQM7/Q13554/Q13555/Q13557
Calculated MW	45,60 kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution::

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.