

Conjugated alpha + beta Synuclein Rabbit mAb

Catalog No: #C56061

Package Size: #C56061-AF350 100ul #C56061-AF405 100ul #C56061-AF488 100ul #C56061-AF555 100ul #C56061-AF594 100ul #C56061-AF647 100ul #C56061-AF680 100ul #C56061-AF750 100ul #C56061-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Conjugated alpha + beta Synuclein Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	alpha + beta Synuclein Antibody detects endogenous levels of total alpha + beta Synuclein
Immunogen Description	A synthesized peptide derived from human alpha + beta Synuclein
Conjugates	AF Dye
Other Names	alpha synuclein; beta synuclein ; NACP; PD1 ; SNCA; SNCB;
Accession No.	Uniprot:P37840/Q16143
Calculated MW	18kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution.;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.