

Conjugated CXCR4 Rabbit mAb

Catalog No: #C56063

Package Size: #C56063-AF350 100ul #C56063-AF405 100ul #C56063-AF488 100ul #C56063-AF555 100ul #C56063-AF594 100ul
#C56063-AF647 100ul #C56063-AF680 100ul #C56063-AF750 100ul #C56063-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Conjugated CXCR4 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse
Specificity	CXCR4 Antibody detects endogenous levels of total CXCR4
Immunogen Description	A synthesized peptide derived from human CXCR4
Conjugates	AF Dye
Other Names	CD184 ; CXCR4; C-X-C chemokine receptor type 4; FB22; Fusin; HM89; LCR1; LESTR; NPYRL; SDF-1 receptor; Stromal cell- derived factor 1 receptor;
Accession No.	Uniprot:P61073
Calculated MW	43kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution;;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.