

Conjugated RUNX1/RUNX2/RUNX3 Rabbit mAb

Catalog No: #C56068

Package Size: #C56068-AF350 100ul #C56068-AF405 100ul #C56068-AF488 100ul #C56068-AF555 100ul #C56068-AF594 100ul
#C56068-AF647 100ul #C56068-AF680 100ul #C56068-AF750 100ul #C56068-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated RUNX1/RUNX2/RUNX3 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	RUNX1/RUNX2/RUNX3 Antibody detects endogenous levels of total RUNX1/RUNX2/RUNX3
Immunogen Description	A synthesized peptide derived from human RUNX1/RUNX2/RUNX3
Conjugates	AF Dye
Other Names	RUNX1; AML1; RUNX2; OSF2; RUNX3; CBFA3;
Accession No.	Uniprot:Q01196/Q13761/Q13950
Calculated MW	Predicted band size: 49 kDa
SDS-PAGE MW	Observed band size: 50 kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution;;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.