

# Conjugated KRAS+HRAS+NRAS Rabbit mAb

Catalog No: #C56091

Package Size: #C56091-AF350 100ul #C56091-AF405 100ul #C56091-AF488 100ul #C56091-AF555 100ul #C56091-AF594 100ul  
 #C56091-AF647 100ul #C56091-AF680 100ul #C56091-AF750 100ul #C56091-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                                                                                                                  |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Conjugated KRAS+HRAS+NRAS Rabbit mAb                                                                                                                                                                             |
| Host Species          | Recombinant Rabbit                                                                                                                                                                                               |
| Clonality             | Monoclonal                                                                                                                                                                                                       |
| Purification          | Affinity-chromatography                                                                                                                                                                                          |
| Species Reactivity    | Human Mouse Rat                                                                                                                                                                                                  |
| Specificity           | KRAS+HRAS+NRAS Antibody detects endogenous levels of total KRAS+HRAS+NRAS                                                                                                                                        |
| Immunogen Description | A synthesized peptide derived from human KRAS+HRAS+NRAS                                                                                                                                                          |
| Conjugates            | AF Dye                                                                                                                                                                                                           |
| Other Names           | GTPase NRas; Transforming protein N-Ras; NRAS; HRAS1; GTPase HRas; H-Ras-1; Ha-Ras; Transforming protein p21; c-H-ras; p21ras; HRAS; HRAS1; GTPase KRas; K-Ras 2; Ki-Ras; c-K-ras; c-Ki-ras; KRAS; KRAS2; RASK2; |
| Accession No.         | Uniprot:P01111/P01112/P01116                                                                                                                                                                                     |
| Calculated MW         | 21kDa                                                                                                                                                                                                            |
| Formulation           | Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.                                                                                                                                                       |
| Storage               | Store at 4ℳ Cin dark for 6 months                                                                                                                                                                                |

## Application Details

Suggested Dilution;;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.