

Conjugated Akt(pan) 1/2/3 Rabbit mAb

Catalog No: #C56093

Package Size: #C56093-AF350 100ul #C56093-AF405 100ul #C56093-AF488 100ul #C56093-AF555 100ul #C56093-AF594 100ul
#C56093-AF647 100ul #C56093-AF680 100ul #C56093-AF750 100ul #C56093-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated Akt(pan) 1/2/3 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Akt(pan)1/2/3 Antibody detects endogenous levels of total AKT1/2/3
Immunogen Description	A synthesized peptide derived from human AKT1/2/3
Conjugates	AF Dye
Other Names	PKB;RAC;PRKBA;PKB-ALPHA;RAC-ALPHA; AKT1 ; AKT2; AKT3; Akt1/2/3; pan Akt;
Accession No.	Uniprot:P31749/P31751/Q9Y243
Calculated MW	56kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution::

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.