

Conjugated IgA Rabbit mAb

Catalog No: #C56140

Package Size: #C56140-AF350 100ul #C56140-AF405 100ul #C56140-AF488 100ul #C56140-AF555 100ul #C56140-AF594 100ul
#C56140-AF647 100ul #C56140-AF680 100ul #C56140-AF750 100ul #C56140-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated IgA Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	IgA Antibody detects endogenous levels of total IgA
Immunogen Description	A synthesized peptide derived from human IgA
Conjugates	AF Dye
Other Names	IGHA1; IGHAI2; PIgR;
Accession No.	Uniprot:P01876/P01877
Calculated MW	60kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution.;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.