

Conjugated Apg10 (Atg10) Rabbit mAb

Catalog No: #C56161

Package Size: #C56161-AF350 100ul #C56161-AF405 100ul #C56161-AF488 100ul #C56161-AF555 100ul #C56161-AF594 100ul
#C56161-AF647 100ul #C56161-AF680 100ul #C56161-AF750 100ul #C56161-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated Apg10 (Atg10) Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	Apg10 (Atg10) Antibody detects endogenous levels of total Apg10 (Atg10)
Immunogen Description	A synthesized peptide derived from human Apg10 (Atg10)
Conjugates	AF Dye
Other Names	APG10L;ATG10 autophagy related 10 homolog (S. cerevisiae); ATG10 autophagy related 10 homolog; Autophagy related protein 10;
Accession No.	Uniprot:Q9H0Y0
Calculated MW	25kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4Λ Cin dark for 6 months

Application Details

Suggested Dilution:;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.