

Conjugated CRISPR-Cas9 SA Rabbit mAb

Catalog No: #C56163

Package Size: #C56163-AF350 100ul #C56163-AF405 100ul #C56163-AF488 100ul #C56163-AF555 100ul #C56163-AF594 100ul
#C56163-AF647 100ul #C56163-AF680 100ul #C56163-AF750 100ul #C56163-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated CRISPR-Cas9 SA Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Staphylococcus aureus
Specificity	CRISPR-Cas9 SA Antibody detects endogenous levels of total CRISPR-Cas9 SA
Immunogen Description	A synthesized peptide derived from CRISPR-Cas9 SA
Conjugates	AF Dye
Other Names	Cas9; CRISPR-associated endonuclease Cas9/Csn1; CRISPR-Cas9/Csn1; csn1; SpyCas9;
Accession No.	Uniprot:J7RUA5
Calculated MW	124kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution.;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.