

Conjugated Drebrin Rabbit mAb

Catalog No: #C56164

Package Size: #C56164-AF350 100ul #C56164-AF405 100ul #C56164-AF488 100ul #C56164-AF555 100ul #C56164-AF594 100ul
#C56164-AF647 100ul #C56164-AF680 100ul #C56164-AF750 100ul #C56164-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated Drebrin Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Drebrin Antibody detects endogenous levels of total Drebrin
Immunogen Description	A synthesized peptide derived from human Drebrin
Conjugates	AF Dye
Other Names	D0S117E; DBN1; Drebrin 1; Drebrin; Drebrin E2; Drebrin1; DrebrinE;
Accession No.	Uniprot:Q16643
Calculated MW	120kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution::

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.