

Conjugated CHD4 Rabbit mAb

Catalog No: #C56169

Package Size: #C56169-AF350 100ul #C56169-AF405 100ul #C56169-AF488 100ul #C56169-AF555 100ul #C56169-AF594 100ul
 #C56169-AF647 100ul #C56169-AF680 100ul #C56169-AF750 100ul #C56169-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	Conjugated CHD4 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	CHD4 Antibody detects endogenous levels of total CHD4
Immunogen Description	A synthesized peptide derived from human CHD4
Conjugates	AF Dye
Other Names	ATP dependent helicase CHD4; CHD4; Chromodomain helicase DNA binding protein 4; Mi 2 autoantigen 218 kDa protein; Mi 2b; Mi2 beta;
Accession No.	Uniprot:Q14839
Calculated MW	Predicted band size: 218 kDa
SDS-PAGE MW	Observed band size: 260 kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution::

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.