

Conjugated Human IgG Rabbit mAb

Catalog No: #C56180

Package Size: #C56180-AF350 100ul #C56180-AF405 100ul #C56180-AF488 100ul #C56180-AF555 100ul #C56180-AF594 100ul
#C56180-AF647 100ul #C56180-AF680 100ul #C56180-AF750 100ul #C56180-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated Human IgG Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	Human IgG Antibody detects endogenous levels of total Human IgG
Immunogen Description	A synthesized peptide derived from human Human IgG
Conjugates	AF Dye
Other Names	G1m marker; G2m marker; G3m marker; G4m marker; HDC; Heavy chain disease protein; Human immunglobulin G; Ig gamma 1 chain C region; Ig gamma 2 chain C region; Ig gamma 3 chain C region; Ig gamma 4 chain C region; IgG heavy chain locus; IGHG1; IGHG2; IGHG3; IGHG4;
Accession No.	Uniprot:P01857/P01859/P01860/P01861
Calculated MW	52kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4Λ Cin dark for 6 months

Application Details

Suggested Dilution;;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.