

Conjugated PRAS40 Rabbit mAb

Catalog No: #C56183

Package Size: #C56183-AF350 100ul #C56183-AF405 100ul #C56183-AF488 100ul #C56183-AF555 100ul #C56183-AF594 100ul
#C56183-AF647 100ul #C56183-AF680 100ul #C56183-AF750 100ul #C56183-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated PRAS40 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse
Specificity	PRAS40 Antibody detects endogenous levels of total PRAS40
Immunogen Description	A synthesized peptide derived from human PRAS40
Conjugates	AF Dye
Other Names	40 kDa proline rich AKT substrate; AKT1 S1; AKT1 substrate 1 (proline rich); AKT1S1; Lobe; PRAS40; Proline rich akt substrate;
Accession No.	Uniprot:Q96B36
Calculated MW	40kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4Λ Cin dark for 6 months

Application Details

Suggested Dilution;;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.