

Conjugated MBD1 Rabbit mAb

Catalog No: #C56200

Package Size: #C56200-AF350 100ul #C56200-AF405 100ul #C56200-AF488 100ul #C56200-AF555 100ul #C56200-AF594 100ul
#C56200-AF647 100ul #C56200-AF680 100ul #C56200-AF750 100ul #C56200-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated MBD1 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	MBD1 Antibody detects endogenous levels of total MBD1
Immunogen Description	A synthesized peptide derived from human MBD1
Conjugates	AF Dye
Other Names	CXXC 3; Mbd1; MECP1 COMPLEX; PCM1; RFT;
Accession No.	Uniprot:Q9UIS9
Calculated MW	67kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution.;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.