

Conjugated AMBP Rabbit mAb

Catalog No: #C56224

Package Size: #C56224-AF350 100ul #C56224-AF405 100ul #C56224-AF488 100ul #C56224-AF555 100ul #C56224-AF594 100ul
#C56224-AF647 100ul #C56224-AF680 100ul #C56224-AF750 100ul #C56224-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated AMBP Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	AMBP Antibody detects endogenous levels of total AMBP
Immunogen Description	A synthesized peptide derived from human AMBP
Conjugates	AF Dye
Other Names	A1M; alpha 1 microglobulin bikunin; Alpha 1 microglycoprotein; AMBP; Bikunin; EDC1; HI30;
Accession No.	Uniprot:P02760
Calculated MW	Predicted band size: 39 kDa
SDS-PAGE MW	Observed band size: 39 kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution.;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.