

## Conjugated CRMP2 Rabbit mAb

Catalog No: #C56227

Package Size: #C56227-AF350 100ul #C56227-AF405 100ul #C56227-AF488 100ul #C56227-AF555 100ul #C56227-AF594 100ul  
#C56227-AF647 100ul #C56227-AF680 100ul #C56227-AF750 100ul #C56227-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                   |
|-----------------------|---------------------------------------------------------------------------------------------------|
| Product Name          | Conjugated CRMP2 Rabbit mAb                                                                       |
| Host Species          | Recombinant Rabbit                                                                                |
| Clonality             | Monoclonal                                                                                        |
| Purification          | Affinity-chromatography                                                                           |
| Species Reactivity    | Human Mouse Rat                                                                                   |
| Specificity           | CRMP2 Antibody detects endogenous levels of total CRMP2                                           |
| Immunogen Description | A synthesized peptide derived from human CRMP2                                                    |
| Conjugates            | AF Dye                                                                                            |
| Other Names           | DRP2; CRMP-2; N2A3; Unc-33-like phosphoprotein 2; ULIP-2; DPYSL2; CRMP2; ULIP2; Musunc33; TOAD64; |
| Accession No.         | Uniprot:Q16555                                                                                    |
| Calculated MW         | 62kDa                                                                                             |
| Formulation           | Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.                                        |
| Storage               | Store at 4A Cin dark for 6 months                                                                 |

## Application Details

Suggested Dilution:;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.