

Conjugated Caspr Rabbit mAb

Catalog No: #C56228

Package Size: #C56228-AF350 100ul #C56228-AF405 100ul #C56228-AF488 100ul #C56228-AF555 100ul #C56228-AF594 100ul
#C56228-AF647 100ul #C56228-AF680 100ul #C56228-AF750 100ul #C56228-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated Caspr Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Caspr Antibody detects endogenous levels of total Caspr
Immunogen Description	A synthesized peptide derived from human Caspr
Conjugates	AF Dye
Other Names	Caspr; Caspr1; CNTNAP; Cntnap1; Contactin associated protein 1; MHDNIV; NCP1; Neurexin-4; Nrnx4; p190; Paranodin;
Accession No.	Uniprot:P78357
Calculated MW	165kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution;;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.