

Conjugated Ran Rabbit mAb

Catalog No: #C56233

Package Size: #C56233-AF350 100ul #C56233-AF405 100ul #C56233-AF488 100ul #C56233-AF555 100ul #C56233-AF594 100ul
#C56233-AF647 100ul #C56233-AF680 100ul #C56233-AF750 100ul #C56233-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated Ran Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Ran Antibody detects endogenous levels of total Ran
Immunogen Description	A synthesized peptide derived from human Ran
Conjugates	AF Dye
Other Names	ARA24; Gsp1; GTPase Ran; Guanosine triphosphatase Ran; LPS; RAN member RAS oncogene family; RanGTPase; Ras like protein TC4; Ras related nuclear protein; Ras-like protein TC4; RASL2 8;
Accession No.	Uniprot:P62826
Calculated MW	24kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution;;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.