

Conjugated Gata6 Rabbit mAb

Catalog No: #C56244

Package Size: #C56244-AF350 100ul #C56244-AF405 100ul #C56244-AF488 100ul #C56244-AF555 100ul #C56244-AF594 100ul
#C56244-AF647 100ul #C56244-AF680 100ul #C56244-AF750 100ul #C56244-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated Gata6 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	Gata6 Antibody detects endogenous levels of total Gata 6
Immunogen Description	A synthesized peptide derived from human Gata 6
Conjugates	AF Dye
Other Names	Gata binding factor 6; GATA-binding factor 6; Gata6; Transcription factor Gata 6;
Accession No.	Uniprot:Q92908
Calculated MW	60kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution::

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.