

Lambda Light chain Conjugated Antibody

Catalog No: #C56257

Package Size: #C56257-AF350 100ul #C56257-AF405 100ul #C56257-AF488 100ul #C56257-AF555 100ul #C56257-AF594 100ul
#C56257-AF647 100ul #C56257-AF680 100ul #C56257-AF750 100ul #C56257-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lambda Light chain Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	Lambda Light chain Antibody detects endogenous levels of total Lambda Light chain
Immunogen Description	A synthesized peptide derived from human Lambda Light chain
Conjugates	AF Dye
Other Names	Bence Jones protein; BJP; Constant region of lambda light chains; Ig lambda chain C regions; IGLC; IGLC1; IGLC2; IGLC3; IGLC6; IGLV; Immunoglobulin lambda constant 1; Mcg marker; Paraprotein;
Accession No.	Uniprot:P01701/P01842
Calculated MW	25kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.