

Conjugated Claudin 4 Rabbit mAb

Catalog No: #C56285

Package Size: #C56285-AF350 100ul #C56285-AF405 100ul #C56285-AF488 100ul #C56285-AF555 100ul #C56285-AF594 100ul
#C56285-AF647 100ul #C56285-AF680 100ul #C56285-AF750 100ul #C56285-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated Claudin 4 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	Claudin 4 Antibody detects endogenous levels of total Claudin 4
Immunogen Description	A synthesized peptide derived from human Claudin 4
Conjugates	AF Dye
Other Names	CLD4; CLDN4; CPE- receptor; CPE-R; CPER; CPETR1; Clostridium perfringens enterotoxin receptor; WBSR8; Williams-Beuren syndrome chromosome region 8 protein;
Accession No.	Uniprot:O14493
Calculated MW	18kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4A Cin dark for 6 months

Application Details

Suggested Dilution;;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.