

TRPM7 Conjugated Antibody

Catalog No: #C56296

Package Size: #C56296-AF350 100ul #C56296-AF405 100ul #C56296-AF488 100ul #C56296-AF555 100ul #C56296-AF594 100ul #C56296-AF647 100ul #C56296-AF680 100ul #C56296-AF750 100ul #C56296-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	TRPM7 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	TRPM7 Antibody detects endogenous levels of total TRPM7
Immunogen Description	A synthesized peptide derived from human TRPM7
Conjugates	AF Dye
Other Names	ALSPDC; CHAK; CHAK1; LTrpC7; TRPM7;
Accession No.	Uniprot:Q96QT4
Calculated MW	213kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

- AF350 conjugated: most applications: 1: 50 - 1: 250
- AF405 conjugated: most applications: 1: 50 - 1: 250
- AF488 conjugated: most applications: 1: 50 - 1: 250
- AF555 conjugated: most applications: 1: 50 - 1: 250
- AF594 conjugated: most applications: 1: 50 - 1: 250
- AF647 conjugated: most applications: 1: 50 - 1: 250
- AF680 conjugated: most applications: 1: 50 - 1: 250
- AF750 conjugated: most applications: 1: 50 - 1: 250
- Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.