

## CD33 Conjugated Antibody

Catalog No: #C56297

Package Size: #C56297-AF350 100ul #C56297-AF405 100ul #C56297-AF488 100ul #C56297-AF555 100ul #C56297-AF594 100ul  
#C56297-AF647 100ul #C56297-AF680 100ul #C56297-AF750 100ul #C56297-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                       |
|-----------------------|-------------------------------------------------------|
| Product Name          | CD33 Conjugated Antibody                              |
| Host Species          | Rabbit                                                |
| Clonality             | Monoclonal                                            |
| Isotype               | Rabbit IgG                                            |
| Purification          | Affinity-chromatography                               |
| Species Reactivity    | Human                                                 |
| Specificity           | CD33 Antibody detects endogenous levels of total CD33 |
| Immunogen Description | A synthesized peptide derived from human CD33         |
| Conjugates            | AF Dye                                                |
| Other Names           | CD33; gp67; My9; p67; SIGLEC3;                        |
| Accession No.         | Uniprot:P20138                                        |
| Calculated MW         | 65-75kDa                                              |
| Storage               | Store at 4°C in dark for 6 months                     |

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.