

ZEB1 Conjugated Antibody

Catalog No: #C56299

Package Size: #C56299-AF350 100ul #C56299-AF405 100ul #C56299-AF488 100ul #C56299-AF555 100ul #C56299-AF594 100ul
#C56299-AF647 100ul #C56299-AF680 100ul #C56299-AF750 100ul #C56299-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZEB1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	AREB6 Antibody detects endogenous levels of total AREB6
Immunogen Description	A synthesized peptide derived from human AREB6
Conjugates	AF Dye
Other Names	BZP; TCF8; AREB6; FECD6; NIL2A; PPCD3; ZFHEP; ZFHX1A; DELTAEF1;
Accession No.	Uniprot:P37275
Calculated MW	200kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.