

# Conjugated GFP Rabbit mAb

Catalog No: #C56313

Package Size: #C56313-AF350 100ul #C56313-AF405 100ul #C56313-AF488 100ul #C56313-AF555 100ul #C56313-AF594 100ul  
 #C56313-AF647 100ul #C56313-AF680 100ul #C56313-AF750 100ul #C56313-Biotin 100ul  
 Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                     |
|-----------------------|-----------------------------------------------------------------------------------------------------|
| Product Name          | Conjugated GFP Rabbit mAb                                                                           |
| Host Species          | Recombinant Rabbit                                                                                  |
| Clonality             | Monoclonal                                                                                          |
| Purification          | Affinity-chromatography                                                                             |
| Species Reactivity    | All                                                                                                 |
| Specificity           | GFP antibody detects endogenous levels of C-terminal, internal, and N-terminal GFP-tagged proteins. |
| Immunogen Description | A synthesized peptide derived from green fluorescent protein                                        |
| Conjugates            | AF Dye                                                                                              |
| Other Names           | GFP ; yfp ; Green fluorescent protein ;                                                             |
| Accession No.         | Uniprot:Q963I9                                                                                      |
| Calculated MW         | 27kDa                                                                                               |
| Formulation           | Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.                                          |
| Storage               | Store at 4A Cin dark for 6 months                                                                   |

## Application Details

Suggested Dilution.;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.