

# Alpha1 adaptin Conjugated Antibody

Catalog No: #C56328

Package Size: #C56328-AF350 100ul #C56328-AF405 100ul #C56328-AF488 100ul #C56328-AF555 100ul #C56328-AF594 100ul  
#C56328-AF647 100ul #C56328-AF680 100ul #C56328-AF750 100ul #C56328-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                           |
|-----------------------|---------------------------------------------------------------------------|
| Product Name          | Alpha1 adaptin Conjugated Antibody                                        |
| Host Species          | Rabbit                                                                    |
| Clonality             | Monoclonal                                                                |
| Isotype               | Rabbit IgG                                                                |
| Purification          | Affinity-chromatography                                                   |
| Species Reactivity    | Human Mouse Rat                                                           |
| Specificity           | Alpha1-adaptin Antibody detects endogenous levels of total Alpha1-adaptin |
| Immunogen Description | A synthesized peptide derived from human Alpha1-adaptin                   |
| Conjugates            | AF Dye                                                                    |
| Other Names           | ADTAA; CLAPA1;                                                            |
| Accession No.         | Uniprot:O95782                                                            |
| Calculated MW         | 108kDa                                                                    |
| Storage               | Store at 4°C in dark for 6 months                                         |

## Application Details

### Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.