

# PEDF Conjugated Antibody

Catalog No: #C56332

Package Size: #C56332-AF350 100ul #C56332-AF405 100ul #C56332-AF488 100ul #C56332-AF555 100ul #C56332-AF594 100ul  
 #C56332-AF647 100ul #C56332-AF680 100ul #C56332-AF750 100ul #C56332-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PEDF Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	PEDF Antibody detects endogenous levels of total PEDF
Immunogen Description	A synthesized peptide derived from human PEDF
Conjugates	AF Dye
Other Names	EPC1; OI12; OI6; PEDF; PIG35; Serpinf1;
Accession No.	Uniprot:P36955
Calculated MW	46kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

### Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.