

IGFBP7 Conjugated Antibody

Catalog No: #C56335

Package Size: #C56335-AF350 100ul #C56335-AF405 100ul #C56335-AF488 100ul #C56335-AF555 100ul #C56335-AF594 100ul
#C56335-AF647 100ul #C56335-AF680 100ul #C56335-AF750 100ul #C56335-Biotin 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	IGFBP7 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	IGFBP7 Antibody detects endogenous levels of total IGFBP7
Immunogen Description	A synthesized peptide derived from human IGFBP7
Conjugates	AF Dye
Other Names	AGM; PSF; TAF; FSTL2; IBP-7; MAC25; IGFBP-7; RAMSVPS; IGFBP-7v; IGFBPRP1;
Accession No.	Uniprot:Q16270
Calculated MW	29kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

- AF350 conjugated: most applications: 1: 50 - 1: 250
- AF405 conjugated: most applications: 1: 50 - 1: 250
- AF488 conjugated: most applications: 1: 50 - 1: 250
- AF555 conjugated: most applications: 1: 50 - 1: 250
- AF594 conjugated: most applications: 1: 50 - 1: 250
- AF647 conjugated: most applications: 1: 50 - 1: 250
- AF680 conjugated: most applications: 1: 50 - 1: 250
- AF750 conjugated: most applications: 1: 50 - 1: 250
- Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.