

# MSMB Conjugated Antibody

Catalog No: #C56355

Package Size: #C56355-AF350 100ul #C56355-AF405 100ul #C56355-AF488 100ul #C56355-AF555 100ul #C56355-AF594 100ul  
#C56355-AF647 100ul #C56355-AF680 100ul #C56355-AF750 100ul #C56355-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                           |
|-----------------------|-----------------------------------------------------------------------------------------------------------|
| Product Name          | MSMB Conjugated Antibody                                                                                  |
| Host Species          | Rabbit                                                                                                    |
| Clonality             | Monoclonal                                                                                                |
| Isotype               | Rabbit IgG                                                                                                |
| Purification          | Affinity-chromatography                                                                                   |
| Species Reactivity    | Human                                                                                                     |
| Specificity           | Prostate Secretory Protein/PSP Antibody detects endogenous levels of total Prostate Secretory Protein/PSP |
| Immunogen Description | A synthesized peptide derived from human Prostate Secretory Protein/PSP                                   |
| Conjugates            | AF Dye                                                                                                    |
| Other Names           | MSP; PSP; IGBF; MSPB; PN44; PRPS; HPC13; PSP57; PSP94; PSP-94;                                            |
| Accession No.         | Uniprot:P08118                                                                                            |
| Calculated MW         | 15kDa                                                                                                     |
| Storage               | Store at 4°C in dark for 6 months                                                                         |

## Application Details

### Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.