

PAR6 Conjugated Antibody

Catalog No: #C56373

Package Size: #C56373-AF350 100ul #C56373-AF405 100ul #C56373-AF488 100ul #C56373-AF555 100ul #C56373-AF594 100ul
#C56373-AF647 100ul #C56373-AF680 100ul #C56373-AF750 100ul #C56373-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PAR6 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	PAR6 Antibody detects endogenous levels of total PAR6
Immunogen Description	A synthesized peptide derived from human PAR6
Conjugates	AF Dye
Other Names	PAR6; PAR6C; TAX40; PAR-6A; TIP-40; PAR6alpha;
Accession No.	Uniprot:Q9NPB6
Calculated MW	37kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.