

Prostatic Acid Phosphatase Conjugated Antibody

Catalog No: #C56384

Package Size: #C56384-AF350 100ul #C56384-AF405 100ul #C56384-AF488 100ul #C56384-AF555 100ul #C56384-AF594 100ul
#C56384-AF647 100ul #C56384-AF680 100ul #C56384-AF750 100ul #C56384-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prostatic Acid Phosphatase Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	PAP (prostatic acid phosphatase) Antibody detects endogenous levels of total PAP (prostatic acid phosphatase)
Immunogen Description	A synthesized peptide derived from human PAP (prostatic acid phosphatase)
Conjugates	AF Dye
Other Names	5'-NT; 5'-nucleotidase; ACP3; PAP; Prostatic acid phosphatase; TMPase;
Accession No.	Uniprot:P15309
Calculated MW	50kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.