

Pancreatic Lipase $\alpha\Omega^{1/2}\alpha\Omega^{1/2}$ PNLIPo $\Omega^{1/2}\alpha\Omega^{1/2}$
Conjugated Antibody

Catalog No: #C56391

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #C56391-AF350 100ul #C56391-AF405 100ul #C56391-AF488 100ul #C56391-AF555 100ul #C56391-AF594 100ul
#C56391-AF647 100ul #C56391-AF680 100ul #C56391-AF750 100ul #C56391-Biotin 100ul

Description

Product Name	Pancreatic Lipase $\alpha\Omega^{1/2}\alpha\Omega^{1/2}$ PNLIPo $\Omega^{1/2}\alpha\Omega^{1/2}$ Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	Pancreatic Lipase $\alpha\Omega^{1/2}\alpha\Omega^{1/2}$ PNLIPo $\Omega^{1/2}\alpha\Omega^{1/2}$ Antibody detects endogenous levels of total Pancreatic Lipase $\alpha\Omega^{1/2}\alpha\Omega^{1/2}$ PNLIPo $\Omega^{1/2}\alpha\Omega^{1/2}$
Immunogen Description	A synthesized peptide derived from human Pancreatic Lipase $\alpha\Omega^{1/2}\alpha\Omega^{1/2}$ PNLIPo $\Omega^{1/2}\alpha\Omega^{1/2}$
Conjugates	AF Dye
Other Names	Pancreatic lipase; PL; PNLIP; PNLIPD; PTL;
Accession No.	Uniprot:P16233
Calculated MW	49kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

- AF350 conjugated: most applications: 1: 50 - 1: 250
- AF405 conjugated: most applications: 1: 50 - 1: 250
- AF488 conjugated: most applications: 1: 50 - 1: 250
- AF555 conjugated: most applications: 1: 50 - 1: 250
- AF594 conjugated: most applications: 1: 50 - 1: 250
- AF647 conjugated: most applications: 1: 50 - 1: 250
- AF680 conjugated: most applications: 1: 50 - 1: 250
- AF750 conjugated: most applications: 1: 50 - 1: 250
- Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.