

Conjugated Placental lactogen Rabbit mAb

Catalog No: #C56394

Package Size: #C56394-AF350 100ul #C56394-AF405 100ul #C56394-AF488 100ul #C56394-AF555 100ul #C56394-AF594 100ul
#C56394-AF647 100ul #C56394-AF680 100ul #C56394-AF750 100ul #C56394-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated Placental lactogen Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human
Specificity	Placental lactogen Antibody detects endogenous levels of total Placental lactogen
Immunogen Description	A synthesized peptide derived from human Placental lactogen
Conjugates	AF Dye
Other Names	CS2; CSA; CSB; CSH1; CSH2; CSMT; hCS; HCSB; Lactogen; mplII; PL; PI2; Placental lactogen; Placental lactogen II; Prl3b1; Prolactin 3B1; rPLII;
Accession No.	Uniprot:P01243/P0DML2/P0DML3
Calculated MW	25kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4Λ Cin dark for 6 months

Application Details

Suggested Dilution;;

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.