

MEF2C Conjugated Antibody

Catalog No: #C56411

Package Size: #C56411-AF350 100ul #C56411-AF405 100ul #C56411-AF488 100ul #C56411-AF555 100ul #C56411-AF594 100ul
#C56411-AF647 100ul #C56411-AF680 100ul #C56411-AF750 100ul #C56411-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MEF2C Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Rat
Specificity	MEF2C Antibody detects endogenous levels of total MEF2C
Immunogen Description	A synthesized peptide derived from human MEF2C
Conjugates	AF Dye
Other Names	C5DELq14.3; DEL5q14.3; MEF2C;
Accession No.	Uniprot:Q06413
Calculated MW	50-60kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.