

Liver Carboxylesterase 1 Conjugated Antibody

Catalog No: #C56413

Package Size: #C56413-AF350 100ul #C56413-AF405 100ul #C56413-AF488 100ul #C56413-AF555 100ul #C56413-AF594 100ul
#C56413-AF647 100ul #C56413-AF680 100ul #C56413-AF750 100ul #C56413-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Liver Carboxylesterase 1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Liver Carboxylesterase 1 Antibody detects endogenous levels of total Liver Carboxylesterase 1
Immunogen Description	A synthesized peptide derived from human Liver Carboxylesterase 1
Conjugates	AF Dye
Other Names	ACAT; CE 1; CEH; CES1; CES2; CESDD1; Egasyn; ES-HTEL; ES-x; Es22; Esterase 22; hCE 1; HMSE; HMSE1; REH; SES1; TGH; Triacylglycerol hydrolase;
Accession No.	Uniprot:P23141
Calculated MW	55kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.