

## eIF4A1 Conjugated Antibody

Catalog No: #C56446

Package Size: #C56446-AF350 100ul #C56446-AF405 100ul #C56446-AF488 100ul #C56446-AF555 100ul #C56446-AF594 100ul  
#C56446-AF647 100ul #C56446-AF680 100ul #C56446-AF750 100ul #C56446-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                           |
|-----------------------|-----------------------------------------------------------|
| Product Name          | eIF4A1 Conjugated Antibody                                |
| Host Species          | Rabbit                                                    |
| Clonality             | Monoclonal                                                |
| Isotype               | Rabbit IgG                                                |
| Purification          | Affinity-chromatography                                   |
| Species Reactivity    | Human Mouse Rat                                           |
| Specificity           | eIF4A1 Antibody detects endogenous levels of total eIF4A1 |
| Immunogen Description | A synthesized peptide derived from human eIF4A1           |
| Conjugates            | AF Dye                                                    |
| Other Names           | DDX2; DDX2A; eIF 4A I; EIF4A; eIF4A I; eIF4A-I; EIF4A1;   |
| Accession No.         | Uniprot:P60842                                            |
| Calculated MW         | 46kDa                                                     |
| Storage               | Store at 4°C in dark for 6 months                         |

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.