

# eIF2B epsilon Conjugated Antibody

Catalog No: #C56450

Package Size: #C56450-AF350 100ul #C56450-AF405 100ul #C56450-AF488 100ul #C56450-AF555 100ul #C56450-AF594 100ul  
#C56450-AF647 100ul #C56450-AF680 100ul #C56450-AF750 100ul #C56450-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                         |
|-----------------------|-------------------------------------------------------------------------|
| Product Name          | eIF2B epsilon Conjugated Antibody                                       |
| Host Species          | Rabbit                                                                  |
| Clonality             | Monoclonal                                                              |
| Isotype               | Rabbit IgG                                                              |
| Purification          | Affinity-chromatography                                                 |
| Species Reactivity    | Human                                                                   |
| Specificity           | eIF2B epsilon Antibody detects endogenous levels of total eIF2B epsilon |
| Immunogen Description | A synthesized peptide derived from human eIF2B epsilon                  |
| Conjugates            | AF Dye                                                                  |
| Other Names           | CACH; CLE; EIF 2B; EIF2B; EIF2B5; EIF2BE; EIF2Bepsilon;                 |
| Accession No.         | Uniprot:Q13144                                                          |
| Calculated MW         | 80kDa                                                                   |
| Storage               | Store at 4°C in dark for 6 months                                       |

## Application Details

### Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.