

P70 S6 Kinase alpha Conjugated Antibody

Catalog No: #C56451

Package Size: #C56451-AF350 100ul #C56451-AF405 100ul #C56451-AF488 100ul #C56451-AF555 100ul #C56451-AF594 100ul
#C56451-AF647 100ul #C56451-AF680 100ul #C56451-AF750 100ul #C56451-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	P70 S6 Kinase alpha Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	P70 S6 Kinase alpha Antibody detects endogenous levels of total P70 S6 Kinase alpha
Immunogen Description	A synthesized peptide derived from human P70 S6 Kinase alpha
Conjugates	AF Dye
Other Names	p70 alpha; P70 beta 1; p70 S6K; P70S6K; P70S6K1; p70S6Kb; PS6K; RPS6KB1; S6K; S6K1; STK14A;
Accession No.	Uniprot:P23443
Calculated MW	70kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.