

Prostatic Acid Phosphatase Conjugated Antibody

Catalog No: #C56459

Package Size: #C56459-AF350 100ul #C56459-AF405 100ul #C56459-AF488 100ul #C56459-AF555 100ul #C56459-AF594 100ul
#C56459-AF647 100ul #C56459-AF680 100ul #C56459-AF750 100ul #C56459-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prostatic Acid Phosphatase Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse
Specificity	Prostatic Acid Phosphatase Antibody detects endogenous levels of total Prostatic Acid Phosphatase
Immunogen Description	A synthesized peptide derived from human Prostatic Acid Phosphatase
Conjugates	AF Dye
Other Names	ACP3; PAP; TMPase;
Accession No.	Uniprot:P15309
Calculated MW	52kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.