

Cardiac Troponin C Conjugated Antibody

Catalog No: #C56464

Package Size: #C56464-AF350 100ul #C56464-AF405 100ul #C56464-AF488 100ul #C56464-AF555 100ul #C56464-AF594 100ul
#C56464-AF647 100ul #C56464-AF680 100ul #C56464-AF750 100ul #C56464-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cardiac Troponin C Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Cardiac Troponin C Antibody detects endogenous levels of total Cardiac Troponin C
Immunogen Description	A synthesized peptide derived from human Cardiac Troponin C
Conjugates	AF Dye
Other Names	tnnc1a; TNC; Tnnc1; TNNI3; Troponin C;
Accession No.	Uniprot:P63316
Calculated MW	18kDa
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.